

Set recommended by DELMAR, M.D.

- ① 50230 BA **HOPKINS® II Forward-Oblique Telescope 30°**, enlarged view, diameter 4 mm, length 18 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red
- 50200 ES **Optical Dissector**, with distal spatula, fenestrated, large, sharp, for use with HOPKINS® II telescopes 30°
- ② 50200 EL **Sheath**, with rotating stopcock, with distal beak 11 mm, for restraining tissue, for use with HOPKINS® II telescopes 30°
- ③ 50410 FB **TAKE-APART® Raspatory**, malleable, spoon-shaped spatula, slightly curved, width 12 mm, working length 15 cm, color code: red-yellow
- ④ 50410 CB **TAKE-APART® Raspatory**, malleable, spatula 12 x 12 mm, working length 15 cm, color code: red-yellow
- ⑤ 50410 EB **TAKE-APART® Raspatory**, malleable, spatula slightly curved, width 8 mm, length 15 cm
- ⑥ 50411 C **TAKE-APART® Elevator**, straight, spatula 12 x 12 mm, length 15 cm
- ⑦ 50411 A **TAKE-APART® Elevator**, straight, width 9 mm, working length 15 cm, color code: black
- ⑧ 50410 AB **TAKE-APART® Raspatory**, malleable, width 9 mm, working length 15 cm, color code: red-yellow
- ⑨ 50411 BG **TAKE-APART® Elevator**, curved, width 6 mm, working length 15 cm
- ⑩ 50312 AG **TAKE-APART® Nerve Hook**, curved, hook angled to left, working length 15 cm, color code: black-yellow
- ⑪ 50312 BG **Same**, hook curved to right
- ⑫ 50230 F **CLICK'liné Dissecting and Grasping Forceps**, dismantling, non-rotating, with 4 locking positions, single action jaws, straight sheath, size 3 mm, length 18 cm
- ⑬ 50244 GW **CLICK'liné Dissecting and Grasping Forceps**, dismantling, non-rotating, with 4 locking positions, double action jaws, robust, horizontal opening, curved sheath, size 5 mm, length 18 cm
- ⑭ 50231 GW **CLICK'liné Scissors**, dismantling, non-rotating, with 4 locking positions, double action jaws, curved, horizontal opening, curved sheath, size 3 mm, length 18 cm
- ⑮ 50239 G **CLICK'liné Scissors**, dismantling, non-rotating, with 4 locking positions, double action jaws, curved, serrated, straight sheath, size 5 mm, length 18 cm
- ⑯ 50203 DR **Coagulation Suction Tube**, L-shaped, hook curved to right, curved sheath, diameter 5 mm, working length 15 cm
- 30804 **Handle with Trumpet Valve**, for suction or irrigation, **autoclavable**, for use with 5 mm coagulating suction tubes and 5 mm suction and irrigation tubes