

Set recommended by GRÖNER, M.D.

- ① 50230 BA **HOPKINS® II Forward-Oblique Telescope 30°**, enlarged view, diameter 4 mm, length 18 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red
- 50200 EK **Telescope Sheath**, with rotating stopcock, with distal beak 5 mm, for restraining tissue, for use with HOPKINS® II telescopes 30°
- ② 50200 EM **Sheath**, with rotating stopcock, with distal beak 8 mm, for restraining tissue, for use with HOPKINS® II telescopes 30°
- ③ 50410 CB **TAKE-APART® Raspatory**, malleable, spatula 12 x 12 mm, working length 15 cm, color code: red-yellow
- ④ 50410 FB **TAKE-APART® Raspatory**, malleable, spoon-shaped spatula, slightly curved, width 12 mm, working length 15 cm, color code: red-yellow
- ⑤ 50312 AG **TAKE-APART® Nerve Hook**, curved, hook angled to left, working length 15 cm
- ⑥ 50312 BG **Same**, hook angled to right
- ⑦ 50244 GW **CLICK'liné Dissecting and Grasping Forceps**, dismantling, non-rotating, with 4 locking positions, double action jaws, robust, horizontal opening, curved sheath, size 5 mm, length 18 cm
- ⑧ 50243 G **CLICK'liné Scissors**, dismantling, non-rotating, with 4 locking positions, double action jaws, curved, serrated, straight sheath, size 5 mm, length 18 cm
- ⑨ 50203 DR **Coagulation Suction Tube**, L-shaped, hook curved to right, curved sheath, diameter 5 mm, working length 15 cm
- 30804 **Handle with Trumpet Valve**, for suction or irrigation, **autoclavable**, for use with 5 mm coagulating suction tubes and 5 mm suction and irrigation tubes