

Set recommended by Prof. HÖNIG, M.D.

- ① 50200 L **Optical Retractor**, with handle, width of distal end 8 mm, with fixation screw for control of penetration depth of endoscope, with Telescope Sheath 50200 LS, for use with HOPKINS® II telescopes diameter 4 mm, length 18 cm
- ② 50230 BA **HOPKINS® II Forward-Oblique Telescope 30°**, enlarged view, diameter 4 mm, length 18 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red
- ③ 58210 AGA **Raspatory**, curved, width 9 mm, working length 15 cm
- ④ 58210 GGA **Raspatory**, curved, arc-shaped spatula, width 9 mm, working length 15 cm
- ⑤ 58210 LGA **Raspatory**, curved, for lifting and dissecting the orbita rim periosteum, distal end of spatula 90° curved upwards, upper part sharp, width 9 mm, working length 15 cm
- ⑥ 50200 DS **V-Tunnel Drill System**, with handle, with three teeth for secure fixation on the cranium, with rotating drill guide, for use with Drills 50200 DB/DK
- ⑦ 50200 DB **Drill**, with penetration stop, diameter 2.1 mm, length 7.5 cm, for use with V-Tunnel Drill System 50200 DS
- ⑧ 50200 DK **Same**, diameter 1.3 mm
- ⑨ 50203 BH **TAKE-APART® Bipolar Forceps**, with connector pin for bipolar coagulation, with irrigation connector, with fine jaws, width 1.5 mm, size 5 mm, length 15 cm
- ⑩ 50241 OS **CLICK'linē Dissector**, single action jaws, curved jaws, with 4 locking positions, triangular-shaped, blunt, with straight sheath, width of spatula 12 mm, size 7 mm, length 14 cm
- ⑪ 50239 GW **CLICK'linē Scissors**, dismantling, non-rotating, with 4 locking positions, double action jaws, serrated, curved, horizontal opening, curved sheath, size 5 mm, length 18 cm
- ⑫ 50173 AM **Suture Passer**, for endoscopic forehead-lifting, size 2.8 mm, length 17 cm