

**Set recommended by KRAUSE-BERGMANN, M.D.**

- ① 50230 BA **HOPKINS® II Forward-Oblique Telescope 30°**, enlarged view, diameter 4 mm, length 18 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red
- 50200 ES **Optical Dissector**, with distal spatula, fenestrated, large, sharp, for use with HOPKINS® II telescopes 30°
- ② 50310 SK **TAKE-APART® Raspatory**, straight, spatula slightly curved upwards, width 8 mm, working length 10 cm, color code: red
- ③ 50310 AB **TAKE-APART® Raspatory**, malleable, width 9 mm, working length 15 cm, color code: red-yellow
- ④ 50412 AG **TAKE-APART® Nerve Hook**, curved, hook curved to left, working length 15 cm, color code: black-yellow
- ⑤ 50410 FB **TAKE-APART® Raspatory**, malleable, spoon-shaped spatula, slightly curved, width 12 mm, working length 15 cm, color code: red-yellow
- ⑥ 50230 F **CLICK'liné Dissecting and Grasping Forceps**, dismantling, non-rotating, with 4 locking positions, single action jaws, straight sheath, size 3 mm, length 18 cm
- ⑦ 50241 OS **CLICK'liné Dissector**, dismantling, non-rotating, single action jaws, curved jaws, with 4 locking positions, triangular-shaped, blunt, with straight sheath, width of spatula 12 mm, size 7 mm, length 14 cm
- ⑧ 50244 GW **CLICK'liné Dissecting and Grasping Forceps**, dismantling, non-rotating, with 4 locking positions, double action jaws, robust, horizontal opening, curved sheath, size 5 mm, length 18 cm
- ⑨ 50239 G **CLICK'liné Scissors**, dismantling, non-rotating, with 4 locking positions, double action jaws, curved, serrated, straight sheath, size 5 mm, length 18 cm
- ⑩ 50245 GW **CLICK'liné Scissors**, dismantling, non-rotating, with 4 locking positions, double action jaws, serrated, curved, horizontal opening, curved sheath, size 5 mm, length 18 cm