

SIMMEN-BRINER Special Rhinoplasty Set

①	488060	Ala Double Hook , with octagonal handle, with 2 sharp points, width 2 mm, length 16.5 cm
②	2x 477098	Ala Dilator , for preparation and cutting of nostril cartilage, graduated, width 5 x 8 mm, length 9.5 cm
③	506400	AUFRICHT Nasal Retractor , length 16.5 cm
④	488054	AUFRICHT Nasal Retractor , with side protector, length 16.5 cm
⑤	488070	JOSEPH Elevator , slightly curved, width 4 mm, length 17.5 cm
⑥	489084	Angular Scissors , extra delicate, length 9.5 cm
⑦	488048	KILNER Scissors , curved, sharp/sharp, length 13.5 cm
⑧	513700	FOMON Scissors , curved, working length 6.5 cm
⑨	498822	AIACH Cartilage Graft Forceps , length 15 cm
⑩	516115	RYDER Needle Holder , extra delicate, length 15 cm
⑪	487012	RUBIN Osteotome , width of cut 12 mm, length 16.5 cm
⑫	487014	Same , width of cut 14 mm
⑬	487016	Same , width of cut 16 mm
⑭	486102 SB	Osteotome , width 2 mm, length 19 cm
⑮	523834	Nasal Rasp , tungsten carbide, Fig. 3 and 4
⑯	523856	Same , Fig. 5 and 6
⑰	523878	Same , Fig. 7 and 8
⑱	200600	LEMPERT Rongeur , straight, length 15.5 cm
⑲	469200	ROWLAND-MEYER Nasal Hump Forceps , double joint, curved blades, length 20 cm
⑳	488038	RUBIN Septum Morcelizer , straight, length 20 cm
㉑	498820	AIACH Grid , size 90 x 50 mm
㉒	231000	Hone , "ARKANSAS" oil stone, wedge-shaped, size 10 x 4 cm
㉓	523913	Metal Tray , size 10 x 15 cm

Accessories recommended for containers (not illustrated):

39502 Z **Wire Tray for Cleaning, Sterilization and Storage**

39740 A2 **Sterilization Container**, with MicroStop®

39763 A2 **Coding Label**

39762 D2 **Color Tab**, blue

Further accessories for cleaning and sterilization see catalog HYGIENE