

BERGHAUS Instrument Set for Septo-Rhinoplasty

- ① 488002 COTTLE **Nasal Speculum**, extra slender, without set screw, special matt finish, blade length 35 mm, length 13 cm
- ② 488003 **Same**, blade length 55 mm
- ③ 488004 **Same**, blade length 75 mm
- ④ 488018 KILLIAN-CLAUS **Septum Chisel**, with V-shaped cutting edge, special matt finish, width 5 mm, length 16.5 cm
- ⑤ 488010 RUBIN **Osteotome**, flat, straight, double-edged grinding, rounded corners, with finger grip for controlled osteotome guidance, special matt finish, width 10 mm, length 16.5 cm
- ⑥ 488012 **Same**, width 12 mm
- ⑦ 488014 **Same**, width 14 mm
- ⑧ 488022 **Chisel**, with finger button, flat, special matt finish, width 2 mm, length 18 cm
- ⑨ 488023 **Same**, width 3 mm
- ⑩ 488024 **Same**, width 4 mm
- ⑪ 488030 **Marking Instrument**, for nasal tip, with ball and half cup, straight, special matt finish, length 11 cm
- ⑫ 488032 HEYMANN **Nasal Scissors**, medium size, special matt finish, working length 9.5 cm
- ⑬ 488034 CRAIG **Septum Forceps**, straight, special matt finish, working length 9 cm
- ⑭ 488041 BLAKESLEY **Nasal Forceps**, slender, straight, special matt finish, size 1, working length 11 cm
- ⑮ 488038 RUBIN **Septum Morcelizer**, double action, straight, special matt finish, length 20 cm
- ⑯ 488036 BECKER-CAPLAN **Septum Scissors**, double action, serrated, slender, special matt finish, working length 9.5 cm
- ⑰ 488044 JANSEN **Nasal Dressing Forceps**, bayonet-shaped, slender, special matt finish, length 16 cm
- ⑱ 488045 ADSON-BROWN **Forceps**, micro-model, atraumatic, fine side grasping teeth, special matt finish, length 12 cm
- ⑲ 488046 ADSON **Tissue Forceps**, atraumatic, fine side grasping teeth, special matt finish, length 15 cm
- ⑳ 488047 ADSON **Tissue Forceps**, 1 x 2 teeth, length 12 cm
- ㉑ 488049 KILNER **Scissors**, curved, with tungsten carbide inserts, flat end, sharp/sharp, extra delicate, special matt finish, length 13.5 cm
- ㉒ 488050 KILNER **Scissors**, curved, flat end, blunt/blunt, special matt finish, length 15.5 cm
- ㉓ 488054 AUFRICHT **Nasal Retractor**, with side protector, special matt finish, length 16.5 cm
- ㉔ 488060 **Ala Double Hook**, with octagonal handle, with 2 sharp points, strongly curved, special matt finish, width 2 mm, length 16.5 cm
- ㉕ 488056 KILNER **Ala Retractor**, two sharp points with button, special matt finish, width 10 mm, length 8.5 cm
- ㉖ 488058 **Same**, width 13 mm
- ㉗ 488084 FERGUSON **Suction Tube**, with cut-off hole and stylet, LUER, special matt finish, outer diameter 8 Fr., working length 11 cm
- ㉘ 488085 **Same**, outer diameter 10 Fr.
- ㉙ 488086 **Same**, outer diameter 12 Fr.
- ㉚ 488065 BERGHAUS **Ala Guiding Set**
- ㉛ 488070 JOSEPH **Elevator**, slightly curved, special matt finish, width 4 mm, length 17.5 cm
- ㉜ 488071 **Same**, strongly curved, width 3.4 mm
- ㉝ 488090 **Scalpel Handle**, No. 3, special matt finish, length 12.5 cm
- ㉞ 488074 FREER **Elevator**, double-ended, sharp and blunt, special matt finish, length 20 cm
- ㉟ 488078 **Nasal Rasp**, double-ended, coarse, special matt finish, length 21.5 cm
- ㊱ 488076 AUFRICHT **Glabella Rasp**, double-ended, backward cutting, length 20.5 cm
- ㊲ 488080 RYDER **Needle Holder**, tungsten carbide inserts, extra delicate and flat, special matt finish, length 18 cm
- ㊳ 439300 BERGHAUS **Metal Nasal Mirror**
- 174200 COTTLE **Metal Mallet**, length 18 cm (not illustrated)
- 231000 **Hone**, "ARKANSAS" oil stone, wedge-shaped, size 10 x 4 cm (not illustrated)

Accessories recommended for containers (not illustrated):

- 39502 Z **Wire Tray for Cleaning, Sterilization and Storage**
- 39751 A2 **Sterilization Container**, with MicroStop®
- 39763 A2 **Coding Label**
- 39762 D2 **Color Tab**, blue

Further accessories for cleaning and sterilization see catalog HYGIENE