

Basic Set for Rhinoplasty

①	496400	MASING Surgical Handle
	208015	Blade , Fig. 15, non-sterile, package of 100 (not illustrated)
	208011	Same , Fig. 11 (not illustrated)
②	496500	Surgical Handle , for miniature blades, round, length 15.5 cm
	496764	Miniature Blade , Fig. 64, round, sterile, package of 24 (not illustrated)
③	496900	MASING Nasal Knife
④	493000	COTTLE Nasal Knife , rounded edge, length 14 cm
⑤	3x 417013	Cotton Applicators , triangular serrated, length 15 cm
⑥	511210	Scissors , extra delicate, curved, length 10 cm
⑦	513410	COTTLE Scissors , curved, length 10.5 cm
⑧	511414	JOSEPH Scissors , curved, length 14 cm
⑨	513200	WALTER Scissors , angled, length 10 cm
⑩	513700	FOMON Dorsal Scissors , angular, working length 6.5 cm
⑪	484004	COTTLE Chisel , flat, width 4 mm
⑫	484007	Same , width 7 mm
⑬	484009	Same , width 9 mm
⑭	484206	COTTLE Crossbar Osteotome , graduated, straight
⑮	484406	Same , curved
⑯	484800	MASING Chisel , straight
⑰	484700	Same , curved to left
⑱	484600	Same , curved to right
⑲	498400	COTTLE Retractor
⑳	505700	COTTLE Retractor , double-ended
㉑	505800	NEIVERT Retractor , double-ended, made of wire
㉒	505100	COTTLE Retractor , 2 prongs, blunt prong on left
㉓	505000	Same , blunt prong on right
㉔	499101	Hook , one prong, sharp curve, length 16.5 cm
㉕	506400	AUFRICHT Nasal Retractor
㉖	479800	Suction Elevator , with stylet
㉗	479000	MASING Elevator , double-ended
㉘	478100	JOSEPH-MASING Elevator
㉙	478305	McKENTY Elevator , width 5 mm
㉚	479200	COTTLE Elevator , double-ended
㉛	534015	COTTLE Lower Lateral Forceps , bayonet-shaped
㉜	533212	ADSON-BROWN Tissue Forceps , atraumatic, fine side grasping teeth, length 12 cm
㉝	523600	Nasal Rasp , double-ended, fine
㉞	523700	Nasal Rasp , double-ended, coarse
㉟	528800	WALSHAM Septum Straightening Forceps
㊱	174200	COTTLE Metal Mallet
㊲	523000	JOSEPH Saw , bayonet-shaped, right
㊳	523100	Same , left
㊴	403655	COTTLE Nasal Speculum , blade length 55 mm, length 13 cm
㊵	403575	KILLIAN-STRUYCKEN Nasal Speculum , standard model, blade length 75 mm, length 13.5 cm
㊶	403590	Same , blade length 90 mm
㊷	456001	BLAKESLEY Nasal Forceps , size 1
㊸	456003	Same , size 3
㊹	458000	WATSON-WILLIAMS Nasal Forceps , serrated
㊺	523900	COTTLE Bone Crusher , with clip
㊻	517000	MASING Needle Holder , ear forceps shape
㊼	516513	NEIVERT-MASING Needle Holder , length 13 cm
㊽	515017	Needle Holder , extra delicate, length 17 cm
㊾	529207	FRAZIER Suction Tube , 7 Fr.
㊿	529209	Same , 9 Fr.
①	529105	Suction Tube , curved, outer diameter 5 mm
②	534500	COTTLE Columella Clamp , length 11 cm
③	469000	KAZANIJAN Nasal Hump Forceps , length 20 cm
④	468000	JANSEN-MIDDLETON Septum Forceps

Accessories recommended for containers (not illustrated):

- 39502 Z **Wire Tray for Cleaning, Sterilization and Storage**
- 39751 A2 **Sterilization Container**, with MicroStop®
- 39763 A2 **Coding Label**
- 39762 D2 **Color Tab**, blue

Further accessories for cleaning and sterilization see catalog HYGIENE