

Basic Instruments for Wrist Arthroscopy



Basic Instruments for Wrist Arthroscopy



Laying the Foundation for Precise and Efficient Work in Hand Surgery

As an anatomically small structure, the hand places great demands on instrument precision in wrist arthroscopy. KARL STORZ has compiled an arthroscopy set and an accessory set that are optimized to the demands of wrist arthroscopy.

The proven KARL STORZ HOPKINS® arthroscopes ensure brilliant visualization even in tight joint spaces.



The instruments of the SILCUT® PRO line are perfectly designed for use in small joints as well. Thanks to the optimal combination of cutting force and cutting sensation, the force can be uniformly applied to the target tissue.



The extremely flat jaws of the SILCUT® PRO line allow successful work on structures even in tight joint areas.

The basic instruments we have compiled for you offer optimal conditions for precise wrist arthroscopy.

Wrist Arthroscopy Set



28320	Wrist Arthroscopy Set consisting of:
28300 BA	HOPKINS® Wide Angle Forward-Oblique Telescope 30° , diameter 2.4 mm, length 10 cm, autoclavable , fiber optic light transmission incorporated, color code: red
28303 DS	High-Flow Arthroscope Sheath , with snap-in coupling mechanism, diameter 3.5 mm, working length 6.5 cm, two stopcocks, rotatable, for use with HOPKINS® Telescope 30° 28300 BA and Obturators 28302 BU/BV
28302 BV	Obturator , blunt, for use with Arthroscope Sheaths 28303 DS/BS
28305 BA	HOPKINS® Forward-Oblique Telescope 30° , diameter 1.9 mm, length 6.5 cm, autoclavable , fiber optic light transmission incorporated, color code: red
28306 BN	High-Flow Arthroscope Sheath , with standard coupling mechanism, diameter 2.5 mm, working length 4 cm, one stopcock, fixed, for use with HOPKINS® Telescope 30° 28305 BA and Obturators 28306 BS/BT, color code: red
28306 BT	Obturator , blunt, for use with Arthroscope Sheaths 28306 BN
28145 SN	Probe , graduated, length of hook 2 mm, diameter 1 mm, working length 4 cm
28302 S	Probe , graduated, length of hook 1 mm, diameter 1 mm, working length 7.5 cm
28572 NJ	SILGRASP® PRO Suture Forceps , straight jaws, sheath diameter 2.8 mm, straight, handle with cleaning connector, working length 8.5 cm
28572 AD	SILCUT® PRO Punch , through-cutting, cutting width 2 mm, straight jaws, sheath diameter 2.5 mm, straight, handle with cleaning connector, working length 8.5 cm
28572 O	SILCUT® PRO Oval Biter , through-cutting, cross-toothed, straight jaws, sheath diameter 2.5 mm, straight, handle with cleaning connector, working length 8.5 cm
28572 ADL	SILCUT® PRO Punch , through-cutting, cutting width 2 mm, jaws curved to left 45°, sheath diameter 2.5 mm, straight, handle with cleaning connector, working length 8.5 cm
28572 ADR	SILCUT® PRO Punch , through-cutting, cutting width 2 mm, jaws curved to right 45°, sheath diameter 2.5 mm, straight, handle with cleaning connector, working length 8.5 cm
221150	HARTMANN Ear Forceps , extra delicate, serrated, 1 x 4.5 mm, working length 12.5 cm
517000	MASING Needle Holder , ear forceps shape, intranasal, smooth jaws, working length 8 cm
215010	Retractor , 2 x 3 prongs, blunt, length 10 cm

SilGrasp® PRO Grasping Forceps



NEW

Sheath diameter 2.8 mm, straight sheath, working length 8.5 cm



28572 CG

SILGRASP® PRO Cartilage Grasping Forceps, spoon-shaped jaws, straight jaws



28572 CGU

SILGRASP® PRO Cartilage Grasping Forceps, spoon-shaped jaws, jaws curved, 30° upwards



28572 GU

SILGRASP® PRO Grasping Forceps, spoon-shaped jaws, serrated, jaws curved, 30° upwards



28572 G

SILGRASP® PRO Grasping Forceps, spoon-shaped jaws, serrated, straight jaws



28572 AG

SILGRASP® PRO Alligator Grasping Forceps, serrated, straight jaws

DRILLCUT-X® ARTHRO Shaver Handpiece



28200 DX

DRILLCUT-X® ARTHRO Shaver Handpiece, up to 8000 rpm, for use with UNIDRIVE® S III ARTHRO SCB, as of software status 1.10.

The shaver handpieces should NOT be treated with universal spray!

Burrs and Shaver Blades, for single use



28206 FAS **Small Joint Burr**, for single use, sterile, diameter 2.5 mm, working length 70 mm, package of 6, for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL SCB, color code: orange



28206 FBS **Small Joint Burr**, for single use, sterile, diameter 3.5 mm, working length 70 mm, package of 6, for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL SCB, color code: orange



28206 DAS **Aggressive Full Radius Resector**, for single use, sterile, diameter 2.5 mm, working length 70 mm, package of 6, for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL SCB, color code: blue



28206 DBS **Aggressive Full Radius Resector**, for single use, sterile, diameter 3.5 mm, working length 70 mm, package of 6, for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL SCB, color code: blue

Reusable Shaver Blades

28206 DA **Aggressive Full Radius Resector**, reusable, sterilizable, diameter 2.5 mm, working length 70 mm, for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL SCB, color code: blue

28206 DB **Aggressive Full Radius Resector**, reusable, sterilizable, diameter 3.5 mm, working length 70 mm, for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL SCB, color code: blue

Wrist Arthroscopy Accessory Set



28330 Wrist Arthroscopy Accessory Set

consisting of:

- 791618 **Ligature Scissors**, with tungsten carbide blades, curved, lower blade serrated, length 18 cm, color code: one black handle ring, one gold-plated handle ring
- 792005 **MAYO Scissors**, curved, length 17 cm
- 791803 **JAMESON Scissors**, curved, pointed, delicate tips, length 15.5 cm
- 200400 **BEYER Rongeur**, slender, curved, length 17 cm
- 516015 **Needle Holder**, tungsten carbide inserts, length 15 cm
- 2x 499003 **Hook**, one prong, large curve, length 18 cm
- 2x 499207 **JOSEPH Double Hook**, sharp, width 7 mm, length 15 cm
- 2x 535312 **HALSTEAD "Mosquito" Artery Forceps**, micro-model, curved, length 12.5 cm
- 535112 **HALSTEAD "Mosquito" Artery Forceps**, micro-model, straight, length 12.5 cm
- 474006 **Double Elevator**, size 18 cm
- 2x 505400 **Retractor**, double-ended, size 6 x 20 mm and 5 x 10 mm, length 14.5 cm
- 533116 **ADSON Tissue Forceps**, 1 x 2 teeth, length 15 cm
- 533022 **ADSON Dressing Forceps**, serrated, tungsten carbide inserts, length 12 cm
- 2x 208000 **Surgical Handle**, Fig. 3, length 12.5 cm, for Blades 208010 – 15, 208210 – 15
- 2x 796513 **Towel Forceps**, with serrated inserts, length 13 cm, for paper cover sheets
- 810810 **Medicine Cup**, 100 cm³, diameter 80 mm, height 30 mm

Note: All products can be ordered separately

Notes

It is recommended to check the suitability of the product for the intended procedure prior to use.



KARL STORZ GmbH & Co. KG
Mittelstraße 8, 78532 Tuttlingen, Germany
Postbox 230, 78503 Tuttlingen, Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com

