

STEINER Basic LASER Laryngoscopy Set for Adults

- ① 8666 AN **STEINER Distending Operating Laryngoscope**, for adults, lightweight construction, channel removing vapor integrated in blade and handle, length 18 cm
- ② 8666 DN **STEINER Distending Operating Oropharyngoscope**, for adults, with right and left side wings to prevent the tongue/soft tissue from obstructing the lumen, channel removing vapor integrated in handle and blade, length 14 cm, for use with Clip 497 AC
- ③ 8661 CN **STEINER Operating Laryngoscope**, medium-large, for adults, with integrated channel to remove vapor, length 18 cm
- ④ 8661 DN **Same**, small, for adults with difficult anatomical circumstances and for children, length 19 cm
- ⑤ 8663 BH **Grasping Forceps**, curved to right, serrated, insulated sheath, working length 23 cm
- ⑥ 8663 CH **Same**, curved to left
- ⑦ 8593 GM **Miniature Grasping Forceps**, extra delicate, serrated, with triangular jaws, curved upwards to right, working length 23 cm
- ⑧ 8593 HM **Same**, curved upwards to left
- ⑨ 8662 D **LARYNGOFORCE®II Grasping Forceps**, delicate, serrated, special matt finish, working length 22 cm
- ⑩ 8662 E **LARYNGOFORCE®II Grasping Forceps**, with oval alligator jaws, small, working length 22 cm
- ⑪ 8662 FL **Same**, medium, working length 22 cm
- ⑫ 8662 GL **LARYNGOFORCE®II Grasping Forceps**, with oval, fenestrated alligator jaws, large, with suction channel to remove LASER vapors, working length 22 cm
- ⑬ 8665 T **Clips**, titanium LT 200, medium, 5 mm, sterile, package of 36 cartridges with 6 clips each
- ⑭ 8665 L **LARYNGOFORCE®II Clip Forceps**, jaw curved to left, working length 22 cm
- ⑮ 8665 R **Same**, jaw curved to right
- ⑯ 8606 D **Coagulation Suction Tube**, with ergonomic handle, insulated, with connection for unipolar coagulation, outer diameter 2 mm, working length 23 cm
- ⑰ 8606 E **Same**, outer diameter 2.5 mm
- ⑱ 8606 F **Same**, outer diameter 3.5 mm
- 26005 M **Unipolar High Frequency Cord**, with 5 mm plug, length 300 cm, for KARL STORZ and Erbe units (not illustrated)
- ⑲ 8596 R **Protector**, to protect tissue against unintentional LASER irradiation, oval, diameter 7 mm, working length 23 cm
- ⑳ 497 AC **Clip**, for proximal illumination
- ㉑ 8575 K **Laryngoscope Holder and Chest Support**, GÖTTINGEN model, rod with metal ring, diameter 9 cm, length 34 cm

If required:

- 8575 L **Support Table**, GÖTTINGEN model, for Laryngoscope Holder 8575 K/KC, height adjustable swivel arm with movable plate, can be mounted on operation table equipped with standard sliding rail 25 x 10 mm (not illustrated)

Accessories recommended for containers (not illustrated):

- 39502 Z **Wire Tray for Cleaning, Sterilization and Storage**
- 39752 A2 **Sterilization Container**, with MicroStop®
- 39763 A2 **Coding Label**
- 39762 D2 **Color Tab**, blue

Further accessories for cleaning and sterilization see catalog HYGIENE