

Basic Panendoscopy Set

- ① 8575 C RIECKER **Laryngoscope**, slotted, adult size, with Fiber Optic Light Carrier 8575 CF
- ② 8590 B KLEINSASSER **Operating Laryngoscope**, large, length 17 cm
- ③ 8590 C KLEINSASSER **Operating Laryngoscope**, medium, length 17 cm (most commonly used model)
- ④ 8590 DN KLEINSASSER **Operating Laryngoscope**, small, length 18 cm (for particularly difficult anatomical circumstances)
- ⑤ 8575 K **Laryngoscope Holder and Chest Support**, GÖTTINGEN model
including:
Laryngoscope Holder, GÖTTINGEN model, with adjustment wheel
Support Rod, movable, with metal ring, diameter 9 cm, length 34 cm
- ⑥ 12060 AV ROBERTS-JESBERG **Esophagoscope**, oval, size 12 x 16, length 40 cm, for use with Fiber Optic Light Carrier 12061 AV and/or Prismatic Light Deflector 10101 FA
- ⑦ 12061 AV **Fiber Optic Light Carrier**, for use with esophagoscopes, length 40 cm
- ⑧ 8590 GF **Fiber Optic Light Carrier**, for distal illumination, working length 14 cm
- 495 UD **Fiber Optic Light Cable**, with straight connector, 2x diameter 3.5 mm, length 180 cm, for simultaneous light transmission for 2 instruments (not illustrated)
- ⑨ 12070 **Handle**
- ⑩ 810820 **Cup Medicine**, 200 ccm, diameter 100 mm, height 50 mm
- ⑪ 810810 **Cup Medicine**, 100 ccm, diameter 80 mm, height 30 mm
- ⑫ 818500 **Medicine Dish**, glass, diameter 32 mm, colored
- ⑬ 744600 McIVOR **Mouth Gag**, with removable bow for replacement of rubber tube
- ⑭ 744303 **Tongue Blade**, WÜRZBURG model, slotted, size 3
- ⑮ 744304 **Same**, size 4
- ⑯ 400500 HARTMANN **Nasal Speculum**, for adults, length 13 cm
- ⑰ 754719 SCHNIDT **Forceps**, slightly curved, serrated, length 19 cm
- ⑱ 426516 JANSEN **Nasal Dressing Forceps**, bayonet-shaped, length 16.5 cm
- ⑲ 204810 FERGUSON **Suction Tube**, with cut-off hole and stylet, LUER, 10 Fr., working length 11 cm
- ⑳ 769905 **Laryngeal Mirror**, with unscrewable hexagonal handle, size 5
- ㉑ 8596 W **Cotton Applicator**, straight, working length 25 cm
- ㉒ 8575 QS **Teeth Protector**, silicone, autoclavable
- ㉓ 8591 A KLEINSASSER **Cutting Forceps**, straight, with 2 mm cupped jaws, working length 23 cm
- ㉔ 8591 B **Same**, curved upwards
- ㉕ 8591 C **Same**, curved to right
- ㉖ 8591 D **Same**, curved to left
- ㉗ 8594 A KLEINSASSER **Scissors**, straight, working length 23 cm
- ㉘ 8594 B **Same**, angled 45°
- ㉙ 8594 C **Same**, curved to right
- ㉚ 8594 D **Same**, curved to left
- ㉛ 8591 P **Cutting Forceps**, with round 5 mm cupped jaws, robust model, double action jaws, sheath diameter 2.5 mm, working length 23 cm
- ㉜ 8596 D **Laryngeal Probe**, straight, sheath angled, round handle, working length 23 cm
- ㉝ 10380 AK **Suction Tube**, diameter 4 mm, working length 35 cm
- ㉞ 10381 B **Suction Tube**, with rubber tip, straight, diameter 2 mm, working length 35 cm
- ㉟ 10381 A **Suction Tube**, with rubber tip, straight, diameter 4 mm, working length 50 cm

Accessories recommended for containers (not illustrated):

- 2x 39502 Z **Wire Tray for Cleaning, Sterilization and Storage**
- 39753 A2 **Sterilization Container**, with MicroStop®
- 39763 A2 **Coding Label**
- 39762 D2 **Color Tab**, blue

Further accessories for cleaning and sterilization see catalog HYGIENE