

Basic Pediatric Laryngoscopy Set

- ① 8574 J BENJAMIN **Operating Laryngoscope**, for children, length 15 cm
- ② 8574 JH **Fiber Optic Light Carrier**, for distal illumination, length 12 cm, for use with 8574 J/JP
- ③ 8587 KK LINDHOLM **Operating Laryngoscope**, for children, distal end width 16 mm, length 11 cm
- ④ 8587 GF **Fiber Optic Light Carrier**, for distal illumination, length 7.5 cm, for use with Laryngoscope 8587 A/AA/KK
- ⑤ 8587 P BENJAMIN-LINDHOLM **Operating Laryngoscope**, for newborn babies and infants, inner diameter of proximal end 26 x 16 mm, distal end 14 x 11 mm, length 9.5 cm
- ⑥ 8587 PF **Fiber Optic Light Carrier**, for distal illumination, length 6.5 cm, for use with 8587 P
- ⑦ 497 AC **Clip**, for proximal illumination
- ⑧ 8602 KS KLEINSASSER **Suction Tube**, ball-shaped end, outer diameter 2 mm, working length 18 cm
- ⑨ 8603 KS **Same**, outer diameter 3 mm
- ⑩ 8591 MK **Cutting Forceps, with cupped jaws**, jaws open horizontally, distal end malleable, diameter 4 mm, working length 18 cm
- ⑪ 8593 AK KLEINSASSER **Grasping Forceps**, straight, serrated, working length 18 cm
- ⑫ 8591 BK KLEINSASSER **Cutting Forceps**, upturned, with cupped jaws, 2 mm, working length 18 cm
- ⑬ 8591 AK **Same**, straight
- ⑭ 8594 DK KLEINSASSER **Scissors**, curved to left, working length 18 cm
- ⑮ 8594 CK **Same**, curved to right
- ⑯ 8594 AK KLEINSASSER **Scissors**, straight, working length 18 cm
- ⑰ 8575 KC BENJAMIN-PARSONS **Laryngoscope Holder and Chest Support**, rod with metal ring, with 2 lateral screws

Accessories recommended for containers (not illustrated):

39502 Z **Wire Tray for Cleaning, Sterilization and Storage**

39752 A2 **Sterilization Container**, with MicroStop®

39763 A2 **Coding Label**

39762 D2 **Color Tab**, blue

Further accessories for cleaning and sterilization see catalog HYGIENE