

# CASTELNUOVO Special Set for Endonasal Endoscopic Anterior Skull Base Surgery

- ① 7230 AA **HOPKINS® II Straight Forward Telescope 0**, diameter 4 mm, length 18 cm
- ② 7230 AS **Irrigation Sheath**, for 7230 AA, for use with CLEARVISION® II System (see chapter 4, SINOSCOPY)
- ③ 7230 FA **HOPKINS® II Forward-Oblique Telescope 45°**, diameter 4 mm, length 18 cm
- ④ 7230 FS **Irrigation Sheath**, for 7230 FA, for use with CLEARVISION® II System (see chapter 4, SINOSCOPY)
- ⑤ 452831 **RHINOFORCE® II Miniature Nasal Forceps**, with extra fine flat jaws, through-cutting, tissue-sparing, straight jaws, width of cut 1.5 mm, working length 13 cm
- ⑥ 452832 **Same**, jaws 45° upturned
- ⑦ 452833 **Same**, sheath curved 30°, straight jaws
- ⑧ 452834 **Same**, jaws 45° angled upwards
- ⑨ 452841 R CASTELNUOVO **RHINOFORCE® II Miniature Nasal Forceps**, with extra fine flat jaws, through-cutting, tissue-sparing, 65° upturned, right side opening, width of cut 1.5 mm, working length 13 cm
- ⑩ 452841 **Same**, backward opening
- ⑪ 452841 L **Same**, left side opening
- ⑫ 451510 B CASTELNUOVO **RHINOFORCE® II Nasal Forceps**, end of sheath 25° upturned, through-cutting, jaws 45° upturned, BLAKESLEY shape, width 3 mm, working length 13 cm
- ⑬ 456511 B BLAKESLEY-CASTELNUOVO **RHINOFORCE® II Nasal Forceps**, end of sheath 25° upturned, jaws 45° upturned, width 3.5 mm, working length 13 cm
- ⑭ 456010 B BLAKESLEY-CASTELNUOVO **RHINOFORCE® II Nasal Forceps**, end of sheath 25° upturned, with straight jaws, width 3 mm, working length 13 cm
- ⑮ 451010 B CASTELNUOVO **RHINOFORCE® II Nasal Forceps**, end of sheath 25° upturned, through-cutting, with straight jaws, BLAKESLEY shape, width 3 mm, working length 13 cm
- ⑯ 456510 B BLAKESLEY-CASTELNUOVO **RHINOFORCE® II Nasal Forceps**, end of sheath 25° upturned, jaw 45° upturned, width 3 mm, working length 13 cm
- ⑰ 653000 HEUWIESER **Antrum Grasping Forceps**, jaws curved downwards, fixed jaw curved 90°, movable jaw opening up to 120°, backward opening, working length 10 cm
- ⑱ 648523 McKENTY-CASTELNUOVO **Sphenoid Punch**, tip curved 30°, through-cutting, reversible, size 1.6 x 2 mm, working length 17 cm
- ⑲ 629824 CASTELNUOVO **Frontal Sinus Probe, double-ended**, curved, length 22 cm
- ⑳ 28164 EC CASTELNUOVO **Elevator, double-ended**, blunt end angled, graduated, length 26 cm
- ㉑ 28164 EB CASTELNUOVO **Elevator, double-ended**, blunt end slightly curved, length 26 cm
- ㉒ 28164 EA CASTELNUOVO **Elevator, double-ended**, semisharp and blunt, length 26 cm
- ㉓ 28164 H CASTELNUOVO **Hook**, 90°, blunt, with round handle, length 25 cm
- ㉔ 28164 MP **Round Knife**, vertical, oval, with round handle, 3.5 x 2.5 mm, length 25 cm
- ㉕ 28164 ES CASTELNUOVO **Suction Elevator**, curved, length 25 cm
- ㉖ 586146 v. EICKEN-CASTELNUOVO **Antrum Cannula**, S-shaped, strongly curved, outer diameter 2.5 mm
- ㉗ 586145 **Same**, slightly curved
- ㉘ 649180 N FERGUSON-CASTELNUOVO **Suction Tube**, diameter 2 mm, working length 15 cm

#### Accessories recommended for containers (not illustrated):

- 2x 39502 Z **Wire Tray for Cleaning, Sterilization and Storage**
- 39751 A2 **Sterilization Container**, with MicroStop®
- 39763 A2 **Coding Label**
- 39762 D2 **Color Tab**, blue

Further accessories for cleaning and sterilization see catalog HYGIENE