

LENARZ Supplement Set for Cochlear Implantation

①	208000	Surgical Handle , Fig. 3, length 12.5 cm
②	220211	PLESTER Retractor , 2 x 2 prongs, length 11 cm
③	219511	WULLSTEIN Retractor , 3 x 3 prongs, length 11 cm
④	2x 214920	ANDERSON-ADSON Retractor , 4 x 4 prongs, sharp, curved, length 20 cm
⑤	662477	OLIVECRONA Brain Spatula , angled on both sides, concave, elastic, width 11 and 13 mm, length 18 cm
⑥	662476	Same , width 7 and 9 mm
⑦	204405	PLESTER Suction Tube , with grip plate, cut-off hole and stylet, LUER, 5 Fr., length 20 cm
⑧	204407	Same , 7 Fr.
⑨	204409	Same , 9 Fr.
⑩	2x 204005	Suction Tube , angular, LUER-Lock, outer diameter 0.5 mm, working length 6 cm
⑪	2x 204007	Same , outer diameter 0.7 mm
⑫	204004	WULLSTEIN Suction Handle , for Suction Tubes 204005 – 204025, 204305 – 204330, with cut-off hole, LUER cone, length 10 cm
⑬	233110	Octagonal Handle , with LUER cone, straight, length 17 cm
⑭	262550	Diamond Burr , barrel-shaped, diameter 5 mm, length 10 mm, for creating a flat bone edge during cochlea implantation surgery, blunt tip for protection of the dura, with smooth shaft diameter 2.35 mm, length 44 mm
⑮	262506	Diamond Burr , shaft diameter 2.35 mm, diameter 0.6 mm, length 70 mm
⑯	262510	Same , diameter 1 mm
⑰	262514	Same , diameter 1.4 mm
⑱	223891	Seeker , extra delicate, angled 25°, with ball end diameter 1 mm, length 16 cm
⑲	223892	Same , with ball end diameter 1.5 mm
⑳	233212	Forceps , for insertion of cochlea electrodes, nucleus 24 contour, curved to left, working length 6 cm
㉑	233211	Same , curved to right
㉒	233105	Micro Guiding Instrument , for electrodes array, claw-shaped, tip angled downwards, length 16 cm
㉓	233106	Same , tip angled upwards
㉔	224802	Hook , footplate, 0.2 mm, length 16 cm
㉕	224806	PLESTER Hook , footplate, 0.6 mm, length 16 cm
㉖	801814	LENARZ Retractor , 43 x 11 mm, slim and flat blade with integrated suction tube, length 21.5 cm
㉗	786500	SCHÖNBORN Thymus Retractor , width 5 mm, length 20.5 cm
㉘	799990	RANEY Applying and Removing Forceps , for scalp haemostasis clips, length 16 cm
	799995	RANEY Scalp Haemostasis Clips , package of 12 (not illustrated)
㉙	516015	Needle Holder , tungsten carbide inserts, length 15 cm
㉚	536909	Dissecting and Ligature Forceps , straight, smooth jaws, length 9.5 cm
㉛	474000	FREER Elevator , double-ended, semisharp and blunt, length 20 cm
㉜	213008	PLESTER Elevator , width 8 mm, length 18 cm
㉝	662721	ALLIS Forceps , with fine teeth, length 20 cm
㉞	214000	WULLSTEIN Forceps , serrated, length 15 cm
㉟	792403	Dressing Forceps , narrow width, length 14.5 cm
㊱	533212	ADSON-BROWN Tissue Forceps , atraumatic, fine side grasping teeth, length 12 cm
㊲	844319	Bipolar Coagulating Forceps , insulated, bayonet-shaped, slender, tip 0.5 mm wide, length 19 cm, for use with Bipolar High Frequency Cords 847000 or 847000 A/E/M/V
	847000 E	Bipolar High Frequency Cord , for KARL STORZ bipolar instruments, for use with AUTOCON® 50/200/350 and Erbe coagulator (not illustrated)
㊳	231008	Preparation Plate , glass, 10 x 8 cm

Accessories recommended for containers (not illustrated):

39502 Z **Wire Tray for Cleaning, Sterilization and Storage**

39740 A2 **Sterilization Container**, with MicroStop®

39763 A2 **Coding Label**

39762 D2 **Color Tab**, blue

Further accessories for cleaning and sterilization see catalog HYGIENE